

DMS 300-4A Universal Holter ECG Recorder



- Multi-Function Holter ECG Recorder
 - 3-Lead Holter ECG (with 4-electrode cable)
 - 12-Lead Holter ECG (with 10-electrode cable)
 - Includes SAECG Late Potential File
 - Transmits Phone ECG (when 54a is released)
- Performs Resting 12-Lead ECG into:
 - Holter PC computer, or any other PC
 - PDA
 - Windows Smart Phone
- Mini-Holter size and weight
- 4,096 Sample Rate
- Multi-Day: 2 to 4 days
- Single AAA battery
- View ECG quality with PC, PDA, or Windows Smart Phone
- No LCD to avoid massive service problems
- Compatible with CardioScan 12
- Width = 50mm; Depth = 14mm; Length = 83mm.
- Weight = 2oz (57g).
- 2-year Warranty